## Coronary Intravascular Lithotripsy (IVL) System

English .....	2
BG Bulgarian/Български .....	5
CS Czech/Čeština .....	9
DA Danish/Dansk .....	12
DE German/Deutsch .....	15
EL Greek/Ελληνικά .....	19
ES Spanish/Español .....	23
ET Estonian/Eesti .....	27
FI Finnish/Suomi .....	31
FR French/Français .....	34
HR Croatian/Hrvatski .....	38
HU Hungarian/Magyar .....	41
IS Icelandic/Íslenska .....	45
IT Italian/Italiano .....	48
LT Lithuanian/Lietuvių k. ....	52
LV Latvian/Latviešu .....	55
NL Dutch/Nederlands .....	58
NO Norwegian/Norsk .....	62
PL Polish/Polski .....	65
PT Portuguese/Português .....	69
RO Romanian/Română .....	73
SK Slovak/Slovenčina .....	76
SL Slovenian/Slovenščina .....	79
SV Swedish/Svenska .....	82

## Coronary Intravascular Lithotripsy (IVL) System Instructions for Use (IFU)

### For use with the Shockwave Medical, Inc. IVL Generator and Connector Cable

#### Indication for Use

The Shockwave C<sup>2</sup> Coronary IVL System is indicated for lithotripsy-enhanced, low-pressure balloon dilatation of calcified, stenotic de novo coronary arteries prior to stenting.

#### Contents: Shockwave C<sup>2</sup> Coronary IVL Catheter (1)

- The following balloon diameters are available: 2.5mm, 2.75mm, 3.0mm, 3.25mm, 3.5mm, 3.75mm, 4.0mm
- Folded Balloon Diameters are: 0.043" max. for 2.5mm to 2.75mm, 0.044" max. for 3.0mm to 3.5mm and 0.046" max. for 3.75mm to 4.0mm
- 12mm balloon working length
- 138cm catheter working length
- 6F guide catheter compatible
- 0.014" (0.36mm) guidewire compatible, (Rx 190cm - 300cm)

#### How Supplied

The IVL Catheter is supplied sterile via e-beam sterilization and is intended for single use only. Do not re-sterilize as this could damage the device and lead to patient injury. Do not reuse the device as this could result in cross-contamination that could result in patient injury. Carefully inspect all packaging for damage or defects prior to use. Do not use the device if there is any sign of breach of the sterile barrier, as this could indicate loss of sterility that could result in patient injury. Do not use the device if there is damage to the package, as this could lead to device malfunction and result in patient injury. Store the IVL Catheter in a cool, dark, dry place. Storage of the device in extreme conditions may damage the device and/or affect device performance that could lead to patient injury.

#### Device Description

The IVL Catheter is a proprietary lithotripsy device delivered through the coronary arterial system of the heart to the site of an otherwise difficult to treat calcified stenosis, including calcified stenosis that are anticipated to exhibit resistance to full balloon dilation or subsequent uniform coronary stent expansion. Energizing the lithotripsy device will generate intermittent sound waves within the target treatment site, disrupting calcium within the lesion and allowing subsequent dilation of a coronary artery stenosis using low balloon pressure. The IVL Catheter comprises integrated lithotripsy emitters for the localized delivery of intermittent sound wave therapy. The system consists of the IVL Catheter, IVL Connector Cable and IVL Generator. The Shockwave C<sup>2</sup> IVL Catheter is available in seven (7) sizes: 2.5x12mm, 2.75x12mm, 3.0x12mm, 3.25x12mm,

3.5x12mm, 3.75x12mm, and 4.0x12mm. The Shockwave C<sup>2</sup> IVL Catheter is compatible with a 6F guiding catheter, has a working length of 138cm, with hydrophilic coating at the distal end, and shaft depth markers at the proximal end. Refer to Figure 1 below for Shockwave C<sup>2</sup> IVL Catheter components.


Figure 1: Shockwave C<sup>2</sup> IVL Catheter

The catheter shaft contains an inflation lumen, a guidewire lumen, and the lithotripsy emitters. The inflation lumen is used for inflation and deflation of the balloon with 50/50 saline/contrast medium. The guidewire lumen enables the use of a 0.014" guidewire to facilitate advancement of the catheter to and through the target stenosis. The system is designed as 'Rapid Exchange' (Rx), so a length (190cm - 300cm) guidewire is indicated. The emitters are positioned along the length of the balloon working length for delivery of lithotripsy therapy. The balloon is located near the distal tip of the catheter. Two radiopaque marker bands within the balloon denote the working length of the balloon to aid in positioning of the balloon during treatment. The balloon is designed to provide an expandable segment of known length and diameter at a specific pressure. The proximal hub has two ports: one for inflation/deflation of the balloon, and one for the connection to the IVL Connector Cable.

#### Required Devices for the Coronary IVL Procedure

The IVL Catheter is to be used exclusively with the IVL Generator, IVL Connector Cable and its accessories. The IVL Connector Cable is a remote actuator which connects the IVL Generator to the IVL Catheter and is used to activate the lithotripsy therapy from the IVL Generator. Refer to the IVL Generator and IVL Connector Cable Operator's Manual for preparation, operation, warnings and precautions, and maintenance of the IVL Generator and IVL Connector Cable.

#### Devices Required But Not Supplied By Shockwave Medical, Inc.

- 6F guide catheter
- 0.014" (0.36mm) Guide Wire (190cm - 300cm Length)
- 5"x96" (13x244 cm) minimum Sterile Sleeve

## Shockwave C<sup>2</sup> Coronary IVL Catheter Balloon Compliance Chart

Pressure	2.5x 12mm	2.75x 12mm	3.0x 12mm	3.25x 12mm	3.5x 12mm	3.75x 12mm	4.0x 12mm
ATM - kPa	Ø (mm)*	Ø (mm)*	Ø (mm)*	Ø (mm)*	Ø (mm)*	Ø (mm)*	Ø (mm)*
4* - 405	2.4	2.6	2.9	3.2	3.4	3.7	3.8
5 - 507	2.4	2.7	2.9	3.2	3.4	3.7	3.8
6* - 608	2.5	2.7	2.9	3.2	3.4	3.7	3.9
7 - 709	2.5	2.7	3.0	3.3	3.5	3.8	3.9
8 - 811	2.5	2.7	3.0	3.3	3.5	3.8	4.0
9 - 912	2.5	2.8	3.1	3.3	3.5	3.9	4.1
10* - 1013	2.5	2.8	3.1	3.4	3.6	3.9	4.1

**Note:** \*Ø (mm) is ± 0.10 mm; 4 ATM is IVL treatment balloon pressure, 6 ATM is nominal balloon pressure and post-treatment pressure and 10 ATM is RBP (Rated Burst Pressure) of the balloon.

#### Shockwave C<sup>2</sup> Coronary IVL System Sequence Chart

The following pulsing sequence must be followed during treatment. Do not utilize a pulsing sequence other than those outlined in the IVL System Sequence Chart below. Insertion of any size Shockwave C<sup>2</sup> IVL Catheter will automatically program the IVL Generator with the following treatment sequence:

Treatment Frequency	1 Pulse per 1 Second
Maximum Number of Continuous Pulses (1 cycle)	10 Pulses
Minimum Pause Time	10 Seconds
Maximum Total Pulses Per Catheter	80 Pulses (8 Cycles)

In the event the user attempts to deliver more than the maximum number of continuous pulses allowed, the IVL Generator is designed to stop automatically. To resume pulsing, wait at least the minimum pause time before resuming therapy. The therapy button must be released and pressed again to resume therapy. For more information, refer to the IVL Generator and IVL Connector Cable Operator's Manual.

The Shockwave C<sup>2</sup> Coronary IVL Catheter will deliver a maximum of 80 pulses or 8 cycles as noted above. If this count is reached, the catheter shall not be used any further. If further therapy is needed, discard this catheter and obtain a new one. **Caution: Do not exceed 80 pulses in the same treatment segment.**

#### Contraindications for Use

The Shockwave C<sup>2</sup> Coronary IVL System is contraindicated for the following:

1. This device is not intended for stent delivery.
2. This device is not intended for use in carotid or cerebrovascular arteries.

## Warnings

1. This device is intended for single (one) time use only. DO NOT re-sterilize and/or reuse.
2. Do not use a device past the expiration date on the label. Use of expired product may result in patient injury.
3. Percutaneous transluminal IVL should only be performed at hospitals with adequate emergency surgical support.
4. Inspect all product components and packaging prior to use. Do not use the device if the device or the packaging has been damaged or if sterility has been compromised. Damaged product could result in patient injury.
5. Do not use the device if the balloon protective sheath cannot be removed easily prior to use. If excessive force is used the catheter could be damaged. Damaged product could result in patient injury.
6. For preparation, operation, warnings and precautions, and maintenance of the IVL Generator and its accessories refer to the IVL Generator Operator's Manual.
7. IVL Connector Cable is non-sterile and must be enclosed in a sterile sleeve prior to and during use.
8. Use only an appropriately sized balloon for the vessel to be treated.
9. Inflate the balloon according to the balloon compliance chart. Balloon pressure should not exceed the rated burst pressure (RBP).
10. Perform all device manipulations under adequate fluoroscopic guidance.
11. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
12. Do not use, or attempt to straighten, a catheter if the shaft has become bent or kinked. Instead, prepare a new catheter.
13. Do not use excessive force or torque on the catheter as this could result in damage to the device components and result in patient injury.
14. Use only one (1) catheter of any given size. If a second catheter of the same size is necessary, DO NOT re-use the first catheter. Discard it before preparing the second catheter.
15. This device should only be used by physicians trained in angiography and intravascular coronary procedures.
16. Physicians must read and understand these instructions prior to use of the device.
17. Treat patients per standard medication or interventional procedures in the event of complications associated with the procedure or device.

## Precautions

1. Use only the recommended 50/50 contrast/saline balloon inflation medium.
2. During the procedure, appropriate anticoagulant therapy must be provided to the patient as needed. Anticoagulant therapy should be continued for a period of time to be determined by the physician after the procedure.

3. Inflate the balloon according to the compliance chart. DO NOT overinflate balloon as this may cause patient injury.
4. Use the IVL Generator in accordance with recommended settings as stated in the IVL Generator Operator's Manual. DO NOT deviate from recommended settings as this may cause patient injury.
5. If an inability to inflate or maintain balloon pressure occurs, remove the catheter and use a new device.
6. If the IVL Catheter appears not to deliver lithotripsy therapy, remove and replace it with another catheter.
7. Ensure that the IVL Catheter is used with a 0.014" guidewire and is inserted through a 6F guiding catheter at least 0.068" (1.72mm) ID. Failure to do so could result in inadequate device performance or patient injury.
8. Precaution should be taken when handling device after exposure to patient, e.g. contact with blood. Used product is considered biohazardous material and should be disposed of properly as per hospital protocol.
9. Precaution should be taken when treating patients with previous stenting within 5mm of target lesion.

## Adverse Effects

Possible adverse effects are consistent with standard catheter-based cardiac interventions and include, but are not limited to, the following:

- Abrupt vessel closure
- Access site pain
- Allergic reaction to contrast medium, anticoagulant and/or antithrombotic therapy
- Arrhythmia
- Arteriovenous fistula
- Bleeding complications
- Cardiac tamponade or pericardial effusion
- Coronary artery occlusion, perforation, rupture or dissection
- Coronary artery spasm
- Death
- Emboli (air, tissue, thrombus or atherosclerotic emboli)
- Emergency or non-emergency coronary artery bypass surgery
- Emergency or non-emergency percutaneous coronary intervention
- Entry site complications
- Fracture of the guide wire or any component of the device that may or may not lead to device embolism, serious injury or surgical intervention
- Hematoma at the vascular access site(s)
- Hemorrhage
- Hypertension/Hypotension
- Infection/sepsis
- Myocardial Infarction
- Myocardial Ischemia or unstable angina
- Peripheral Ischemia
- Pseudoaneurysm
- Renal failure/insufficiency

- Restenosis of the treated coronary artery
- Shock/pulmonary edema
- Stroke
- Vascular complications which may require surgical repair (conversion to open surgery)

In addition, patients may be exposed to other risks associated with coronary interventional procedures, including risks from conscious sedation and local anesthetic, the radiographic contrast agents used during angiography, the drugs given to manage the subject during the procedure, and the radiation exposure from fluoroscopy.

Risks identified as unique to the device and its use:

- Allergic/immunologic reaction to the catheter material(s) or coating
- Device malfunction or failure
- Excess heat at target site
- Premature ventricular contraction

## Procedural Steps

**Caution:** Refer to the IVL Generator and IVL Connector Cable Operator's Manual for preparation, operation, warnings and precautions, and maintenance of the IVL Generator and IVL Connector Cable.

## Preparation

1. Prepare the insertion site using standard sterile technique.
2. Achieve preferred vascular access and place a guidewire and guide catheter.
3. Select a lithotripsy balloon catheter size that is 1:1 based on balloon compliance chart (above) and reference vessel diameter.
4. Remove the IVL Catheter from the package.
5. Prepare the lithotripsy balloon using standard technique. Fill a 20cc syringe with 5cc of 50/50 saline/contrast medium. Attach syringe to inflation port on catheter hub. Pull vacuum at least 3 times, releasing vacuum to allow the fluid to replace the air in the catheter.
6. Fill deflator device with 10cc of 50/50 saline/contrast medium. Disconnect syringe and connect deflator to inflation port of catheter hub ensuring no air is introduced to the system.
7. Remove the protective sheath and shipping mandrel from the IVL Catheter. **Warning:** Do not use the device if the protective sheath or shipping mandrel are difficult to remove or cannot be removed.
8. Flush the guidewire port with saline.
9. Wet the lithotripsy balloon and distal shaft with sterile saline in order to activate the hydrophilic coating.
10. Insert the IVL Connector Cable into the sterile sleeve or probe cover.

- Remove the cap from the proximal end and attach the IVL Catheter Connector (see Fig 1) to the IVL Connector Cable.
- Attach the other end of the same IVL Connector Cable to the IVL Generator.

**Caution:** Care must be taken to avoid applying lithotripsy therapy, i.e. pressing the therapy button of the IVL Connector Cable while lithotripsy balloon is dry and/or un-inflated, as this may damage the balloon.

#### Delivering the Shockwave C<sup>2</sup> IVL Catheter to the Treatment Site

- Using a guidewire and guiding catheter, position catheter proximal of the treatment site.
- Advance the 0.014" guidewire across the treatment site.
- Select a balloon catheter size that is 1:1 with the reference vessel diameter.
- Load the IVL Catheter over the exchange length (190 - 300cm) 0.014" guidewire and through a guiding catheter and advance balloon to the treatment site.
- Position the lithotripsy balloon at the treatment site using the marker bands to aid in positioning.

#### Treating the Site with Intravascular Lithotripsy

- Once the IVL Catheter is in place, record position using fluoroscopy.
- If position is incorrect, adjust the lithotripsy balloon to the correct position.
- Inflate lithotripsy balloon to 4.0 atm.
- Deliver IVL treatment sequence per the Shockwave C<sup>2</sup> IVL Catheter Sequence Chart by pressing the therapy button on the IVL Connector Cable.

- Inflate lithotripsy balloon to reference size per balloon compliance chart and record lesion response on fluoroscopy.
- Following IVL treatment, deflate lithotripsy balloon and wait 30 seconds to re-establish blood flow.
- Repeat steps 3, 4, 5, and 6 to complete a single treatment of 20 pulses.
- Additional treatments can be performed if deemed necessary. If multiple inflations are required due to a lesion length greater than the lithotripsy balloon length, the recommended balloon overlap is at least 2 mm to prevent geographic miss. However, care must be taken not to exceed 80 pulses maximum in the same treatment segment.
- Perform a completion arteriogram to assess post intervention result.
- Deflate the device and confirm that the balloon is fully deflated prior to removing the IVL Catheter.
- Remove the IVL Catheter. If there is difficulty in removing the device through the hemostatic valve due to the lubricity, gently grasp the IVL Catheter with sterile gauze.
- Inspect all components to ensure that the IVL Catheter is intact. If a device malfunction occurs or any defects are noted on the inspection, flush the guidewire lumen and clean the outer surface of the catheter with saline, store the IVL Catheter in a sealed plastic bag, and contact Shockwave Medical, Inc. at [complaints@shockwavemedical.com](mailto:complaints@shockwavemedical.com) for further instructions.

**Note:** IVL generates mechanical pulses and may cause singular ectopic impulses on the electrocardiogram. In patients with implantable pacemakers and defibrillators, these ectopic impulses may interact with the sensing capabilities. Monitoring of the electrocardiographic rhythm and continuous arterial pressure during IVL treatment is required.

**Caution:** IVL Catheter once pulled out of the body should not be reinserted for additional inflation or lithotripsy treatments. Balloon can be damaged in the process.

#### Patient Information

Physicians should instruct patients to seek medical attention immediately for signs and symptoms of recurrent ischemic heart disease. There are no known limitations to normal daily activities. Patients should be instructed to comply with the medication regimen as prescribed by their physician.

#### Return of Devices

If any portion of the Shockwave IVL System fails prior to or during a procedure, discontinue use and contact your local representative and/or email [complaints@shockwavemedical.com](mailto:complaints@shockwavemedical.com).

Patents: [www.shockwavemedical.com/patents](http://www.shockwavemedical.com/patents)

Symbol	Definition
	Do not re-use
	Use by date
	Sterilized using irradiation
	Caution
	Manufacturer
	Do not use if package is damaged
	Keep dry
	Authorized Representative in the European Community
	Keep away from heat
	Batch code
	Catalogue number
	Do not resterilize

Symbol	Definition
	Non-pyrogenic
	Consult instructions for use
	Contains 1 unit (Contents: 1)
	Recommended Guidewire
	Recommended Guide Catheter
	Balloon Diameter
	Balloon Working Length
	Catheter Working Length (Usable Length, UL)
	Conformité Européenne
	Patents. Refer to <a href="http://www.shockwavemedical.com/patents">www.shockwavemedical.com/patents</a>
	Intravascular Lithotripsy

Shockwave Medical, Inc.  
5403 Betsy Ross Drive  
Santa Clara, CA 95054, USA  
[www.shockwavemedical.com](http://www.shockwavemedical.com)

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MedPass SAS