


A VERSATILE CALCIUM-CRACKING TOOL Count On It



COMMON FEMORAL ARTERY DISEASE

Pre-Treatment Angiogram



Diameter Stenosis = 90%
Lesion Length = 50 mm

IVL Treatment




Shockwave M5+: 8.0 mm x 60 mm,
300 pulses

IVL Treatment




Post-IVL



Diameter Stenosis = 10%


BILATERAL ILIAC ARTERY DISEASE

Pre-Procedure Angiogram



Right Common Iliac Diameter Stenosis = 90%
 Left Common Iliac Diameter Stenosis = 85%
 Right External Iliac Diameter Stenosis = 100%

IVL Treatment



Shockwave M5+: 7.0 mm x 60 mm
 Left Common Iliac = 90 pulses

IVL Treatment



Right Common Iliac = 120 pulses
 Right External Iliac = 90 pulses

Post-IVL



Right Common Iliac Diameter Stenosis = 15%
 Left Common Iliac Diameter Stenosis = 10%
 Right External Iliac Diameter Stenosis = 5%

Case courtesy of Stefano Fazzini, MD

SHOCKWAVE M5+ PERIPHERAL IVL CATHETER SPECIFICATIONS

Catalog Number	Diameter (mm)	Balloon Length (mm)	Sheath Compatibility	Catheter Working Length (cm)	Pulses/Cycle	Cycle	Pulses (Max)	Balloon Crossing Profile (in)
M5PIVL7060	7.0	60	6 F *	135	30	10	300	.068
M5PIVL8060	8.0	60	7 F	135	30	10	300	.074

*6F Compatible with Terumo Pinnacle® Destination® Guiding Sheath and Cook Flexor® Ansel Guiding Sheath. Referenced trademarks are trademarks of their respective owners or holders.

Peripheral IVL Important Safety Information

In the United States: Rx only.

Indications for Use—The Shockwave Medical Intravascular Lithotripsy (IVL) System is intended for lithotripsy-enhanced balloon dilatation of lesions, including calcified lesions, in the peripheral vasculature, including the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries. Not for use in the coronary, carotid or cerebral vasculature.

Contraindications—Do not use if unable to pass 0.014" (M5, M5+, S4, E8) or 0.018" (L6) guidewire across the lesion—Not intended for treatment of in-stent restenosis or in coronary, carotid, or cerebrovascular arteries.

Warnings—Only to be used by physicians who are familiar with interventional vascular procedures—Physicians must be trained prior to use of the device—Use the generator in accordance with recommended settings as stated in the Operator's Manual.

Precautions—use only the recommended balloon inflation medium—Appropriate anticoagulant therapy should be administered by the physician—Decision regarding use of distal protection should be made based on physician assessment of treatment lesion morphology.

Adverse effects—Possible adverse effects consistent with standard angioplasty include—Access site complications —Allergy to contrast or blood thinner—Arterial bypass surgery—Bleeding complications—Death—Fracture of guidewire or device—Hypertension/Hypotension—Infection/sepsis—Placement of a stent—renal failure—Shock/pulmonary edema—target vessel stenosis or occlusion—Vascular complications. Risks unique to the device and its use—Allergy to catheter material(s)— Device malfunction or failure—Excess heat at target site.

Prior to use, please reference the Instructions for Use for more information on indications, contraindications, warnings, precautions and adverse events.

www.shockwavemedical.com/IFU

© 2025 Shockwave Medical Inc., All rights reserved. SPL 65346 Rev. F.