

SHOCKWAVE | M⁵

Intravascular Lithotripsy (IVL) for Safe and Effective Treatment of Vascular Calcium



MINIMIZE TRAUMA | **OPTIMIZE** OUTCOMES | **SIMPLIFY** PROCEDURES

Broad Applications in Calcified Peripheral Arteries


Common Femoral Disease

Enable Endovascular CFA Treatment



Pre-Treatment
Angiogram


Images courtesy of Bill Miller, MD



Final Angiogram
(Post IVL + DCB)

Femoropopliteal Disease

Safe and Effective Vessel Prep or Treatment



Pre-Treatment
Angiogram

Images courtesy of Marianne Brodmann, MD



Final Angiogram
(Post IVL)

Symptomatic Iliac Disease

Preservation of Future Treatment Options



Pre-Treatment
Angiogram

Images courtesy of Patricia Furey, MD



Final Angiogram
(Post IVL)

TAVI, TEVAR & EVAR Access

Facilitate Delivery of Large-Bore Devices



Pre-Treatment
Angiogram



Images courtesy of David Caparrelli, MD



Final Angiogram
(Post IVL & 23Fr TEVAR)



IVL GENERATOR AND CONNECTOR CABLE SPECS

Power	110-240 VAC; 50-60Hz; Single Phase, 15A service	 <p>IVL Generator CATALOG NUMBER: IVLGCC</p>
Size	11" (28.0 cm) high x 6" (15.2 cm) wide x 11.5" (29.2 cm) deep	
Weight	15 pounds (6.8 kg)	
Output	Proprietary pulse delivery system. Output voltage 3000 volts peak, pulse frequency 1Hz	
Mobility	Product is designed to be mounted to an IV pole	
Length	5 ft (1.53m)	 <p>IVL Connector Cable CATALOG NUMBER: IVLCC</p>
Compatibility	Proprietary male key distally designed to connect only to catheter.	
Operation	Lithotripsy pulsing is activated by pushing a button on the Connector Cable.	
Use	Re-usable	

M5 CATHETER SPECS

M5 Branding*	Catalog Number Original Branding*	Diameter (MM)	Length (MM)	Sheath Compatibility	Working Length	Pulses/ Cycle	Cycles	Pulses (Max)	Crossing Profile
M5IVL3560	M732LPBC3560DX1	3.5	60	6F	110	30	10	300	.054
M5IVL4060	M732LPBC4060DX1	4.0	60	6F	110	30	10	300	.057
M5IVL4560	M732LPBC4560DX1	4.5	60	6F	110	30	10	300	.058
M5IVL5060	M732LPBC5060DX1	5.0	60	6F	110	30	10	300	.062
M5IVL5560	M732LPBC5560DX1	5.5	60	6F	110	30	10	300	.064
M5IVL6060	M732LPBC6060DX1	6.0	60	6F	110	30	10	300	.066
M5IVL6560	M732LPBC6560DX1	6.5	60	7F	110	30	10	300	.068
M5IVL7060	M732LPBC7060DX1	7.0	60	7F	110	30	10	300	.073

*Products are equivalent; either Catalog Number may be used when ordering and either product may be received.

**Discover how you can treat calcium more effectively with the
Peripheral Intravascular Lithotripsy (IVL) System.**

Visit shockwavemedical.com or call +1877-77-LITHO (+1 877-775-4846) for more information.

Prior to use, please reference the Instructions for Use for more information on indications, contraindications, warnings, precautions and adverse events.

Contact Shockwave Medical at customerservice@shockwavemedical.com.

www.shockwavemedical.com

© 2019 Shockwave Medical Inc. All rights reserved. SPL 62800 Rev. A